

iZoom Headlight

CALL 1-800-763-6901
or visit www.dentlight.com

See Brighter - See Better - See Sharper

Challenging Any Lights in Power and Size

Patent-pending LED zoom beam technology delivering up to 90000 lux cool white light with variable beam sizes.

Variable Focused and Guided Beam

Unlike other headlights, iZoom has no glares with its unique beam quality giving your patient the most comfort.

Versatile and Constant Power Headlamp

With unique heat sink technology and modular headlamp designs, iZoom is cool, bright and easy to service or upgrade.

Smart and Robust Battery

All aluminum light weight (80 grams) and rechargeable battery sets three power levels and prevents cross-contaminations with easy to disinfect battery pouch.

Ease of Use and Adaptation

Fit to headgear, safety shield, flip-up loupes and standard TTL loupes, iZoom headlight is universal in use and adaptation.



Headlamp

Intensity	Level 3	70000 lux
	Level 2	44000 lux
	Level 1	25000 lux
Weight	19 grams	
Diameter	23 mm	
Temperature	6000 K	
Beam size	Variable 1.8 - 3.5 inches	



TTL loupe adapter



Battery Pouch with Belt/Shirt Clip

Battery Pack

Weight	87 grams
Type	Li-Poly
Size	45x76x15
Run Time	

Level 3	4 hours
Level 2	5.5 hours
Level 1	10 hours

Recharging Time 2 hours

Order information

8500010 iZoom headlight system
Includes headlamp, two rechargeable battery packs, power adapter, articulating arms, loupe adapter, and battery pouch.

- 8500020 LED headlamp
- 8500030 Rechargeable battery pack
- 8500040 Power adapter

iZoom. Expect more for less

