

Nano Loupe Light

When size and power are all that matters

The smallest overall light on the market

Focused Beam

Focused with minimum glares
Unique beam quality for the most comfort.

Constant Power

Unique heat sink technology
Constant power throughout the battery run time.

Reliable

With smart non-memory battery, cooling design, flexible wire and plug-start power driver, DentLight Nano works any time you need it.

Smart and Robust Battery

Solid aluminum rechargeable battery pack weighs only 87 grams with seamless backlit membrane buttons that set three power levels. Two battery packs allow continuous use without interruption. Quick and smart built-in charger never overcharges or undercharges and extends the battery life. Adjustable clip-on battery pouch prevents cross-contaminations and allows easy access on belt, shirt pocket, or arm band.

Ease of Use and Adaptation

Fits to all major loupes including Orascope, Zeiss, Designs for Vision, Surgitel and iZoom.

Headlamp

Intensity	Level 3	45000 lux
	Level 2	31000 lux
	Level 1	16000 lux
Weight	6 grams	
Diameter	18 mm	
Temperature	6000 K	
Beam size	2 inches @ 14"	

Battery Pack

Weight	87 grams	
Type	Li-Poly	
Size	45x76x15mm	
Run Time	Level 3	12 hours
	Level 2	7 hours
	Level 1	5 hours
Recharging Time	2 hours	

Order information

8600110 Nano Loupe Light Kit
Including LED headlamp with articulating arms and power cord, two rechargeable battery packs, universal loupe adapters (to prescription glasses, safety glass, and TTL), L-adapter (to Orascope TTL and Zeiss), U-adapter (to Flip-up), mounting screws, clip-on battery pouch, composite filter, filter tool, universal power adapter.

8600111 Nano Economy Kit includes LED light with one rechargeable battery pack, composite filter, mounting adapter of choice (L-adapter or U-adapter to Flip-up), mounting screws, clip-on battery pouch, universal power adapter.



Orascope TTL L-adapter



Zeiss L-adapter



Flip-up U-adapter



TTL Universal Clip-on

1-year warranty
30-day money back guarantee

Copyright 2011 DentLight Inc.

CALL 1-800-763-6901
or visit www.dentlight.com