


iZoom Loupes

SEEING IS BELIEVING

Easy universal inter-pupil distance adjustment
Antireflection coating for minimum glare

iZoom Flip-up loupes are affordable high performance and lightweight magnification devices providing high resolution and high contrast images for dentists, hygienists and dental assistants. The universal loupes make them easy to adapt for general procedures improving accuracy, efficiency and reducing fatigue and chronicle back pains.



Part Numbers

DESCRIPTION

8500050-XXX_XC

Sports Frame

8500051-XXX_XC

Metal Frame

8500052-YYYY_XC

Titanium Variable Magnification

85000053-400

Titanium 4.0X



Magnification XXX: 200 - 2.0X; 250 - 2.5X; 300 - 3.0X

Magnification YYYY: 2530 - 2.5 to 3.0; 3035 - 3.0 to 3.5 variable magnification

Working Distance: _X: S (short), R (regular), L(long), XL (extra long)

		S	R	L	XL
2X; 2.5X; 3X	Working Distance(inch)	11-14	14-18	17-21	20-24
2X	Field of View (inch)	4	4	4	
2.5X	Field of View (inch)	3	3	3	
3X	Field of View (inch)	2.4	2.4	2.4	
3.5X	Field of View (inch)	2	2		
4.0X	Field of View (inch)	2	2		

Frame Color: C: Black, Red, Blue, Silver, Pink

Standard kit comes with box containing iZoom loupes on selected frame, flip paddle, lens covers, fiber cleaning cloth, and screw driver.

CALL 1-800-763-6901
or visit www.dentlight.com

Copyright 2011 DentLight Inc.

DentLight Inc.
1411 E. Campbell Rd. #500,
Richardson, TX 75081, USA
sales@dentlight.com
www.dentlight.com