

# Nano Loupe Light



### Snap-on Filter

Easy flip-down to avoid composite cure

### Mounting Adapter

Easy attaching to existing loupes

### Wire Clip

Secure wire routing

Actual Size

### Lightweight Headlamp

Small penny size  
Focused spot beam  
Heat insulated front barrel

**Small**

**Slim**

**Strong**

### Smallest and Lightest Battery Pack

Easy intensity digital control  
Two battery packs for nonstop run



# DentLight

# Nano

*"Ideal combination of size and power in clinical illumination", Dr. George Freedman*  
*"This has been the best product I've used", DPS evaluator*

## Amazing Cordless Power At Half the Size!

The award winning and clinically proven cordless Nano outputs an amazing beam of 45000 lux focused power with three digital brightness levels. The small size and superb brightness outshine competitive products at just half the price! .

- Cordless with rechargeable smart lithium battery
- Unique focused cool white beam with minimum glares
- Easily removable composite filter
- Up to 12 hours continuous battery run time
- Durable metal case battery construction
  - Easy cleaning and disinfection
- Raised membrane keypad for easy digital power adjustment
- Cordless with direct power plug that eliminates downtime
- Swivel protective battery pouch for easy placement

## Ease of Use and Adaptation

Fit to all major loupes (e.g. Orascope, Zeiss, Designs for Vision, and Surgitel) with five custom loupe adapters

## Incredible Engineering Design

### • Maximum intensity at the smallest size

Level I	16000 lux
Level II	31000 lux
Level III	45000 lux

\* at 12 inch distance

### • Slim and ergonomic

	Weight (grams)	Size (mm)
Headlamp	6	18 (Diameter)
Battery pack	87	45x76x15

### • Flexible wire and wire routing clips

### • Run time:

Level	I	II	III
(hours)	12	7	5

\* 2 battery packs


\* Continuous run with plug-in cord

- Easy to use membrane control buttons
- Three digital presets with audible beep
- Audible warning when battery is low
- **Recharging time:** 2 hours

## Order information

**8600110 Nano Loupe Light Kit**  
LED headlamp with articulating arm, flex power cord, battery packs (2), a choice for one of the five loupe mounting adapters, battery pouch, removable composite filter, and universal power adapter.

**8600111 Nano Economy Light Kit** LED headlamp with articulating arm, flex power cord, battery pack (1), a choice for one of the five loupe mounting adapters, battery pouch, removable composite filter, and universal power adapter.



## Non-memory Li-Poly Battery You Can Count On!

- AC option allows continuous use without interruption
- Smart built-in charger never over/under charges for longer battery life
- Constant output throughout battery run time
- Multiple charge options
  - Direct plug in
  - Space saving Dual-Battery Counter Charging Stand (optional)



DentLight Inc  
1411 E. Campbell Rd. Suite 500  
Richardson, TX 75081  
[www.dentlight.com](http://www.dentlight.com)

**CALL 800-763-6901**  
**or EMAIL [sales@dentlight.com](mailto:sales@dentlight.com)**

**DentLight**  
Setting New Standards Through Innovation