

New Product Release

October 30, 2013

Nano 2S/2C Loupe Light

DentLight Inc., an LED optics innovator in dental device design and manufacturing, is pleased to launch Nano 2C and 2S – **the smallest, lightest and longest running headlights** for professional medical illumination. Building on the achievements of its award-winning predecessor, the Nano 2C and Nano 2S set the pace for a new era of high performance, smart loupe lights.

With patent-pending optics and the most advanced technology design, the Nano 2 series delivers a focused and uniform light beam of 40,000 lux high power for optimized line-of-sight illumination. Small, sleek and smart, the two new models are made of solid aircraft aluminum with drop-proof reliability and are super easy to disinfect.

The Nano 2 series upgrades the power cord with a Teflon cord, which reduces the wire thickness by 30% while improving its mechanical strength ten times. Nano 2C maintains the smallest battery pack size while doubling the battery running time. Nano 2S slightly increases the battery pack thickness while providing the longest battery running time of up to 24 hours.

This small light amazingly displays a battery icon and has an audible beep 1 minute ahead of the low battery shutdown. The streamlined smart digital battery has one large ON/OFF button for easy hands free activation with an elbow.



Nano 2C and 2S are comfortable to fit on any eyewear and are the heavy-duty headlights that dental professionals can count on to save procedure time, improve quality, and efficiency.

Call **800-763-6901** or visit www.dentlight.com