


FUSION⁵ Curing Light

The Most Powerful and Dependable Curing Light Platform!

3 seconds cure and much more ...



FUSION⁵ Bluephase Valo

Focused Dispersed

- 3-second cure with focused beam
- Flexible platform for oral exams



FUSION Curing Light is a 21st century light that is versatile, cost-effective, and valuable in everyday practice. Adding to the value is the ability to conduct oral cancer examination.

Dr. Howard Glazer



The Fusion Light is one GREAT all around light for the dental office. You get the benefits of three lights in one! Curing, Transillumination and Oral Cancer Screening. The Curing power is remarkable.

Dr. John Comisi




Version 5

DentLight

FUSION⁵ Platform

Amazing Power and Style in the Palm of Your Hand!

- Unique focused beam and advanced backlit LCD display:
 - 4000 mW/cm² Plasma Power mode for fast/deep cure
 - 2000 mW/cm² (Pulsed or per light head in Twinhead)
 - Built-in diagnostic and treatment modes
 - Battery and thermal bars for smart control
- Wide spectrum to cure all materials or for photodynamic treatment
 - 420 - 490nm (FUSION⁵)
 - 385 - 490nm (FUSION⁵ Plus/Grand)
- Sleek metal keypad with large activation and program buttons
- Modular design with easy upgrade and service parts
- Cordless with direct power plug that eliminates downtime
- 90 minutes fast recharge time
- Super slim and ergonomic
 - Tip height only 11.5mm
 - Weighs only 3.9oz
 - Measures at 7.8 inch



Built to Last


Robust anodized aluminum metal construction distributes low LED heat, prevent scratches, discoloration and damage from accidental falls.

Smart Lithium-Ion Battery You Can Count On!


- Cordless power with 70 minutes plasma power run time
 - 1400 3-second cures
- Constant output throughout battery run time
- Common platform for soft tissue laser upgrades

Lowest Profile


Easily cure posterior teeth with small 90° head




Curing Light Dispersion



Optical Spectrum



Add-on Functions



PROFIT FROM the Award-winning FUSION platform

Order Information

Full Range Fast Curing Options

PN:7800080 FUSION⁵ Curing Light Kit:

Including Standard Handpiece, Charging Stand, Curing Caps (tacking, 9mm and 13mm), Power Adapter, Light Shield, Barrier Sleeves (100), and Instructions for Use.

PN:7820080 FUSION Plus Curing Light Kit:

Including Plus Handpiece (10mm lens), Charging Stand, Curing Caps (tacking, 9mm and 13mm), Power Adapter, Light Shield, Barrier Sleeves (100), and Instructions for Use.

PN:7830080 FUSION Grand Curing Light Kit:

Including Plus Handpiece (12mm lens), Charging Stand, Power Adapter, Light Shield, Barrier Sleeves (100), and Instructions for Use.

Clinician's Choice, An independent, non-profit, dental education and product testing foundation, Clinicians Report®. Go to March 2019, Volume 12 Issue 3.



Combination Curing and Diagnostic Options

7700512 FUSION-DOE Kit

Including FUSION⁵, Violet Light Head, Green and Orange Fluorescence Goggles, Orange Caries Fluorescence Shield, Curing Caps, Charging Stand, Power Adapter, Barrier Sleeves, and DVD or USB training guide, and Metal Carrying Box.

7710010 FUSION-DOE Deluxe Kit

Including FUSION⁵, Violet and White Light Heads, Green and Orange Fluorescence Goggles, Orange Caries Fluorescence Shield, Curing Caps, Charging Stand, Power Adapter, Barrier Sleeves, (100), and DVD or USB training guide, and Metal Carrying Box.

Add-on Options and Accessories

- 7500152 FUSION⁵-DOE starter Kit
- 8500123 Green Fluorescence Goggle
- 7500088 Orange Fluorescence Goggle
- 7700192 Orange Fluorescence Shield
- 7700092 Green Fluorescence Shield
- 7500060U Violet Light Head
- 7500082 Barrier Sleeves (100/pack)
- 7701100 Green Camera Filter (fit on SLR)
- 7701110 Orange Camera Filter (fit on SLR)
- 8500122 Green Fluorescence Loupe Filter Pairs
- 8501122 Orange Fluorescence Loupe Filter Pairs



A phenomenal company that consistently develops products to meet dentists needs, Dr. Howard Glazer



1825 Summit Ave. #210
Plano TX 75074
Tel: 972-889-8857
Email: sales@dentlight.com
www.dentlight.com