

Lights, Camera, Action!

Wearable Products that Make a Difference Everyday



Advanced Hands Free Platform

The Smartest and Smallest Light

Black Fluorescence Imaging

SafeLoupe Laser Filter

HD Camera

HD Loupes



Nano Freedom is crammed with innovative technologies. It is bright, light, and feels just right. Best of all no stray wires, no inconvenience, no hassle. iOCam provides super sharp magnified live video for enhanced patient communication/referrals, insurance documentation as well as staff and student training/education.

Dr. George Freedman, Richmond Hill, ON, Canada



iZoom Loupes - Very comfortable, lightweight, and cool-looking. Excellent optics, wide field and very good resolution.

Dr. Michael Miller, Houston, TX



The SafeLoupe system is a "must have" product when doing laser procedures.

Dr. Howard Glazer, Fort Lee, NJ

DentLight

Nano Freedom²

Wireless = Effortless



Mounting Adapter
Easy Dovetail Attachment

Lightweight battery
Proximity sensor On/Off
Slide to existing frame temple
Little extra weight on nose or ear

Freedom Lamp
Penny size 5 grams
Focused beam (~3" field)
Easy clip-on filter

Reconfigurable

Nano Technology Teflon Connection
Field Reconfigurable for Function, Fit, and Service

Superior Performance

20000 (one battery) or 40000 lux (two batteries)
3 (one battery) or 6 hours (two batteries)
Color temperature (5500-6500K)

Safeloupe Insert
Fluorescent Filter




Black Lamp
Dual wavelengths

Seeing is Believing!

Nano Black

**Dual Color: White and Violet for
Illumination and Fluorescence**

- Composites
- Caries
- Plaque
- Cancer



Visualization of dental caries and decays

Detecting composites and excess adhesives

Amazing Wireless Light !

NO MORE TANGLED CORDS!

Nano 2S Plus



Freedom Lamp
Super Lightweight

Reconfigurable
Nano Technology Teflon Wire
Field Reconfigurable
Easy Upgrade to Wireless
Memory wire conforming to frame contour

Incredible Engineering Design

- Three Power Levels (40, 20, 16k lux)
- Light Weight (lamp 5g, battery 118g)
- Super Small (lamp 16mm)
- Flexible Thin wire
- 10-24 hours Non-stop run
- Recharging time: 90% in 2 hours

Xtreme²



Surgical precision

Super high contrast
60000 lux

Seeing beyond the surface

Long Lasting Smart Battery

10 to 24 hours running time (3 power settings)
Super reliable MicroUSB connection
Hands free raised elbow activation button
Constant output power throughout run time
Pass drop proof test



DentLight

iOCam HD Video Camera



Custom Fit with
IZoom HD
Nano Black

Custom Surgical Kit
Wireless iOCam
Wireless Xtreme² light

**HD 1280x720
Magnified Images**

Empower you with iOCam
Magnified Imaging and Communication System
Your handsfree assistant
Capture what you see
Share your findings with patients etc.

**Video or
Picture**

Ease of Use and Adaptation

- Fits to loupes (Orascoptic, DFV, Zeiss etc)
- Compatible with most headlights
- Easy slide in and out of the eyewear
- Built-in microphone

Incredible Engineering Design

- **Magnifications at 16-inch distance**
2.0X, 3.0X, or 5.0X
- **Resolution** 1280x720 @ 30 fps
- **Lightweight and small size**

	Weight (grams)	Size (mm)
Imager	4	23 (D)x13 (L)
Lens	8	15 (D)x13 (L)

• Run time:

USB Limited by hard drive
Wifi Battery 4.5 hours (Optional)

- Live recording or still image snapshot
- Compatible with existing patient software including Dexis Imaging and Vixwin



Communicate better
Raise your acceptance rates
Market your practice with an edge

DentLight

Fit on Any Loupes or Eyewears

Add-on Options and Accessories

8600120 Nano Freedom lamp



8600220 Nano Black lamp



8600022 Nano composite filter



8610110 Nano Freedom battery pod



8600530 Nano Freedom battery (2/pk)

8600130-2 Nano 2S Plus Battery Pack

8601140 Nano Freedom Dual Battery Charger

8500224 Female Dovetail adapter



8500124 H-adapter



8500126 Flex-adapter



8501124 T-adapter



8701100 Universal-H combo adapter



8500121 Uc-adapter



8510120 Articulating arm



8510121 Single side arm



8510122 Branch arm



8510123 Camera branch arm



8501032 2S+ MicroUSB cord



8501022 Freedom cord



Dual Battery Charger always keeps a fresh backup battery

DentLight

PROFIT FROM the Award-winning DentLight

HOW TO ORDER

DentLight Headlight

860011 **5** - **2**

Nano Light Type Gen

3=Nano 2S Plus
5=Nano Freedom
6=Nano Freedom Plus
8=Xtreme
9=Nano Black

=Gen1
2=Gen2

Light Kit includes LED lamp, composite filter, battery pod*, battery**, loupe mounting adapter, charger and accessories

Optional metal case.

Lamp 2-year warranty
* 1 (Nano Freedom); 2 (Nano Freedom Plus)
** 1 (Nano 2S Plus); 2 (Nano Freedom, Xtreme);
4 (Nano Freedom Plus, Nano Black)
Note: Nano Black includes SafeLoupe Fluorescent Insert.



DentLight Camera

880011 **0**

iOcam Type

0=Regular
1=Wireless

iOCam kit includes HD camera with USB wire, 3.0x optical lens, fitting adapter, software, and metal case.
Wireless kit additional: Wifi battery, MicroUSB-USB adapter, charger, IOS/Android App.

Warranty: 1 year

DentLight Loupes

850005 **0** - **250** **R** **IP** **C** **A**

Loupe Type Magnification Working Distance Reading IP Distance Case Type Frame

0=Flipup* 250=2.5X S=10-14" 50=Flipup S=Soft Case A=Adilite
2=Variable* 300=3.0X R=14-18" 54=54mm M=Metal Case D=D2
3=Mini* 350=3.5X E=18-22" 56=56mm M=Metal M=Metal
4=TTL 2530=2.5-3.0X 60=60mm R=Rudilite (Ti-temple)
3035=3.0-3.5X 62=62mm Y=Rudy Project Rydon
400=4.0X



Warranty: optics 5 years

DentLight SafeLoupe - Loupe Filter

851 **0** **810**

SafeLoupe Type Wavelengths

0=Safeloupe Clipon
1=Safeloupe Insert

810=810, 940, 980, 1064 (OD4+)
2780=940 2780 (OD4+)
2940=1064 2940 (OD4+)
4500=380-500 (OD2) Orange
4501=380-500 (OD2) Green

A phenomenal company that consistently develops products to meet dentists needs, Dr. Howard Glazer

DentLight