

New Product Release

Nano Freedom² Loupe Light

Plano, Texas – Nov. 13, 2019 - DentLight Inc. announced today the release of Nano Freedom² reconfigurable wireless loupe light that fits on all major loupes and eyewear.

The product is an upgrade of DentLight's wearable headlights. A small LED lamp (5 grams) has two built-in connectors to link to either one battery pod for 3 hours non-stop operation, or two battery pods for 6 hours, or optionally to a battery pack (Nano 2S Plus) for 10 hours operation.

Unlike other wireless headlights that bear the battery weight on the operator's nose making it very uncomfortable over a long working day, Nano Freedom² has the unique battery design that fits its battery on any frame temples. This allows the battery to disperse its weight evenly across the head resulting in much more comfortable user experience. In addition, the battery pod uses proximity sensor for touch-free activation and easy operation. The system is designed for durability, flexibility and minimum recurring cost for heavy duty daily uses.



About DentLight

DentLight has been designing and manufacturing award-winning dental products in Plano, Texas for the past 14 years. Our mission is to assist clinicians in delivering better procedural outcomes through superior optics and optoelectronic devices. We do this by providing innovative product designs, particularly LED lighting based dental equipment including curing lights, oral exam systems (cancer screening and caries detection), headlights, loupes, video cameras, and eye protection kits.

Call 800-763-6901 or visit www.dentlight.com