

FUSION Photobiomodulation Therapy

Wireless is Effortless

Analgesia

Wound Healing

Deep Tissue Repair

Inflammation and Pain Reduction
superficial musculoskeletal injuries

Interchangeable Light Probes



PBM Probe (LLLTT)

Main Control

Battery Assembly

High Power Modular Probe: 1 cm diameter

- Emit dual wavelengths 660nm/810nm: 200mW each at the same time
- Collimated and uniform beam for deep penetration and perfect optical dosage

Amazing Indestructible Design

- High performance sapphire optical lens
- Solid aircraft aluminum for drop proof reliability
- Advanced backlit LCD with precise control over energy dosage
- Two power levels: normal and high
- Sleek metal keypad with large activation and program buttons
- Super slim and ergonomic
 - Tip height only 11.5mm for easy intraoral access
 - Weighs only 3.9oz; Measures at 7.8 inch long

Smart Lithium-Ion Battery

- Cordless power with 100 minutes run time
- Constant output throughout battery run time
- Common platform for 3-sec curing, caries detection, cancer screening, and soft tissue surgical laser
- 90 minutes fast recharge time

PN:7910080 FUSION PBM Light Kit:

Including PBM Handpiece with dual 660/810 wavelengths, Charging Stand, Power Adapter, Light Shield, Barrier Sleeves (100), and Instructions for Use

A phenomenal company that consistently develops products to meet dentists needs, Dr. Howard Glazer

DentLight

1825 Summit Ave. #210
Plano TX 75074
Tel: 972-889-8857
Email: sales@dentlight.com
www.dentlight.com

