

# Nano Black Dual Light



Safeloupe Insert  
Fluorescent Filter

## Black Dual Lamp

Two color LEDs  
Penny size 5 grams  
Focused beam (~3" field)  
Easy clip-on composite filter

## Lightweight battery


Proximity sensor On/Off  
Slide to existing frame temple  
Little extra weight on nose or ear

## Custom Frame Adapter

Fit on any eyewear or loupes

## Dual Color Light: White and Violet

## Hands Free Illumination



Nightvision adhesives

Visualization of  
dental caries and decays

Visual illumination with white  
Enhanced fluorescence with violet

- Adhesives  
- Caries

*Stands out in violet*

## Superior Performance

Kit includes 4 batteries

3 hours non-stop run time

20000 lux super strong light

2 back up batteries recharge in 1 hour

10x stronger fluorescence than handheld black light

*Seeing is Believing!*

Wireless is Effortless

[www.dentlight.com](http://www.dentlight.com)

**DentLight**