

FUSION PBM

Wireless High Power Cold Laser Therapy Light

- Dual Wavelengths 660nm and 810nm
- 400mW or 800mW (Plus) focused beam
- Impact-resistant Solid Aircraft Aluminum
- Collimated and uniform beam for deep penetration and perfect optical dosage



TMD Treatment

Deep Tissue Repair


Promote Wound Healing

Reduction of pain and inflammation

Reduced need for post-op analgesia

Upgrade to FUSION Plus PBM:

- double the power from 400 to 800 mW
- increase head size from 10mm to 12mm with additional 34mm expander for BOTH Intra and Extra oral use



DentLight

1825 Summit Ave. #210
Plano TX 75074
Tel: 972-889-8857
Email: sales@dentlight.com
www.dentlight.com


Award-winning FUSION platform

Powerful and Versatile


Innovative Design

- Ergonomic, intuitive and easy to use
- New sapphire lens with 90° angle for easy intra-oral access
- Seamless integration with backlit intuitive LCD interface
- 360° swivel head for easy access to full oral cavity
- Smart Lithium battery for 60+ minutes constant output
- Modular platform for add-on functions


Add-on Functions
with Interchangeable Head




Twinhead for bulk cure
PN 7810060



FUSION Grand/FUSION Plus
7830060/7820060



Ultrafast Laser
7300100



Caries/Cancer Detection
PN 7500152



Transillumination
7800060W



PBM/PDT
7800060RIR

Order Information

PN:7910080 FUSION PBM Light Kit:

Including FUSION PBM Handpiece, Barrier Sleeves (100), Safety Goggle, Charging Stand, Power Adapter, Clinical Application Guide, and Instructions for Use
2-year warranty (1-year battery)

PN:7920080 FUSION Plus PBM Light Kit:

Including FUSION Plus PBM Handpiece, 34mm Beam Expander, Barrier Sleeves (100), Safety Goggle, Charging Stand, Power Adapter, Clinical Application Guide, and Instructions for Use
2-year warranty (1-year battery)



DentLight

1825 Summit Ave. #210
Plano TX 75074
Tel: 972-889-8857
Email: sales@dentlight.com
www.dentlight.com