


# SafeLoupe

## Adaptive Laser Filter



**Protect Your Eyes  
from Dental Lasers**

**DentLight**

# SafeLoupe™

## Convert your loupes to laser loupes!

### Are you tired and annoyed with clip-in filters?

- Fighting to fit onto your loupes
- Discomfort on your face from
  - Weight
  - Fogging
  - Hurting your eye/face
- Low light and color distortion
- Procedure interruption
- Cross contamination

### Quick and easy adaptation to your loupes

Preconfigured to your loupes in adapter size and the number of articulating arms  
 Facial viewing field laser protection from diode lasers direct and reflected beams.




### Why buy laser loupes if you can convert your own?

#### Try patent-pending SafeLoupe Laser Filter:

Adaptive filter to your loupes

- Clip-on to loupes ocular rim
- Super lightweight and comfortable
- Minimum light loss
- No color distortion
- Hassle free laser procedures
- Custom preconfiguration to your loupes
- Extraocular shield for full protection



### Specifications

- 810 nm OD5
- 940 nm OD5
- 980 nm OD5
- 1064 nm OD4



### Order Information

Part Number	Description
8510810	Safeloupe Diode Laser Filter Including laser filters (2), pre-configured loupe adapter/articulating arms, and extraocular shield (1)
8510450	Safeloupe Curing Light Filter Including curing filters (2), pre-configured loupe adapter/articulating arms.

Please specify your loupe brand/model and laser brand in your purchase.



CALL 800-763-6901  
 or EMAIL [sales@dentlight.com](mailto:sales@dentlight.com)



DentLight Inc.  
 1825 Summit Ave. #210  
 Plano, TX 75074, USA  
[www.dentlight.com](http://www.dentlight.com)