

iOCam HD Video Camera



Custom Fit with
iZoom HD
Nano Black



Custom Surgical Kit
Wireless iOCam
Wireless Xtreme² light

**HD 1280x720
Magnified Images**

Empower you with iOCam
Magnified Imaging and Communication System
Your handsfree assistant
Capture what you see
Share your findings with patients etc.

**Video or
Picture**

Ease of Use and Adaptation

- Fits to loupes (Orascope, DFV, Zeiss etc)
- Compatible with most headlights
- Easy slide in and out of the eyewear
- Built-in microphone

Incredible Engineering Design


- **Magnifications at 16-inch distance**
2.0X, 3.0X, or 5.0X
- **Resolution** 1280x720 @ 30 fps
- **Lightweight and small size**

	Weight (grams)	Size (mm)
Imager	4	23 (D)x13 (L)
Lens	8	15 (D)x13 (L)

• Run time:

USB	Limited by hard drive
Wifi Battery	4.5 hours (Optional)

- **Live recording or still image snapshot**
- **Compatible with existing patient software including**
Dexis Imaging and Vixwin



Communicate better
Raise your acceptance rates
Market your practice with an edge

DentLight