

ULTRAFAST Diode Laser

Wireless is Effortless

Amazing Laser for Everyday Practice

No More Setup Time and Messy Cords

3W Continuous Power

5W Pulsed Power



Make it Ultrafast

Hemostasis

Frenectomy

Gingivectomy

Aphthous Ulcers

Incision & Excision

Root Canal Therapy

Periodontal Sulcular Debridement

Gingival Troughing for Impression

DentLight

1825 Summit Ave. #210

Plano, TX 75074, US

Tel: 800-763-6901

Fax: 972-346-6550

Email: sales@dentlight.com

Ultrafast Diode Laser

Amazing Power at Your Fingertips!

Sleek Wireless Design

- Patent-pending **810nm** wireless diode laser:
 - up to **3W** CW Power for ultrafast soft tissue sculpting
 - up to **5W** Pulsed Power for precision cut
 - Battery and thermal bars for smart control
- Easy activation with large metal keypads
- Minimal office footprint
 - Ergonomic and lightweight at only **3.5oz**
 - Super slim at 6.4 inches lengths

Enhanced Workflow

The compact design and intuitive LCD interface makes laser procedures quick and effortless, boosting productivity and encouraging frequent use

Built to Last

Aircraft aluminum durable construction distributes heat and prevents scratches, discoloration, or damage from accidental drops

Unique Laser Shield

Adjustable Eye Protection for Operators line of sight

Smart Lithium Battery

- 45+ minutes run time at nominal power
- 90 minute fast recharge time
- Constant output throughout battery run time
- Multiple charge options:
 - Space saving charging stand
 - Cordless with direct power plug backup
- Field proven for nearly 20 years (compatible with FUSION curing light)

PN: 7300100 Ultrafast Diode Laser Kit

Including Laser Handpiece, Fiber Tips (5), Barrier Sleeves (50)
Laser Goggles (2), Laser Shield, Charging Stand, and Power Adapter
2-year warranty (1-year battery)

Distinctive Green Aiming Beam

Enhanced visibility of treatment site

Disposable Fiber Tips

Self-bendable
Time Saving
Prevents cross-contamination



Ultrafast Fiber Tips (25/box)

7300192 Ultrafast Regular 5mm Tips | 400µm
7300193 Ultrafast Perio 7mm Tips | 400µm
7300194 Ultrafast Endo 20mm Tips | 200µm

Call 800-763-6901 or Email sales@dentlight.com